



Natural Relief

The Devon51 is a lymphedema compression therapy pump used to help relieve the symptoms associated with chronic lymphatic and venous disorders in the extremities. When the device is powered on, the attached garments use sequential air compression to move lymphatic fluid throughout the body, mimicking the natural flow. This movement can reduce swelling.

Specifications

Weight: 5.5lbs
Inflation Cycle: 60s
Rest Time: 12s
Dimensions: 4.5" x 9.5" x 10"
Cycle Time: 15s/chamber
End Therapy Vacuum Time: 180s
Treatment Time: 0.5 - 3hrs
Pressure Range: 20-80mmHg
Electrical: 120V~,60Hz,Max.45VA
Fuse Rated: 250VAC, 1.0AMP,
SLO-BLO Operation Mode:
Continuous Operation
Chambers: 4
E0651



Safe



Lightweight



Comfortable

Lower Extremity		Model	Length	Circumference		
Size	Type	4 Chamber	Inseam	Ankle	Knee	Thigh/Groin**
Small	Full Leg	D-300S	23in. / 58.4cm.	21in. / 53.3cm.	25in. / 63.5cm.	30in. / 76.2cm.
Small Wide	Full Leg	D-300SW	26in. / 66cm.	24in. / 61cm.	29in. / 73.7cm.	35in. / 88.9cm.
Small Extra Wide	Full Leg	D-300SXW	24in. / 61cm.	28in. / 71.1cm.	35in. / 88.9cm.	40in. / 101.6cm.
Medium	Full Leg	D-300M	29in. / 73.7cm.	22in. / 55.9cm.	27in. / 68.6cm.	31in. / 78.8cm.
Medium Wide	Full Leg	D-300MW	29in. / 73.7cm.	25in. / 63.5cm.	31in. / 78.7cm.	35in. / 88.9cm.
Medium Extra Wide	Full Leg	D-300MXW	28in. / 71.1cm.	29in. / 73.7cm.	36in. / 91.4cm.	40in. / 101.6cm.
Large	Full Leg	D-300L	34in. / 86.4cm.	21in. / 53.3cm.	27in. / 68.6cm.	31in. / 78.7cm.
Large Wide	Full Leg	D-300LW	36in. / 91.4cm.	22in. / 55.9cm.	31in. / 78.7cm.	40in. / 101.6cm.
Half	Half Leg	D-301H	16.5in. / 41.9cm.	22in. / 55.9cm.	28in. / 71.1cm.	–
Half Wide	Half Leg	D-301HW	19in. / 48.3cm.	24in. / 61cm.	30in. / 76.2cm.	–

Upper Extremity		Model	Length*	Circumference		
Size	Type	4 Chamber	Axilla-Fingertips	Axilla	Elbow	Wrist
Medium	Arm	D-302M	26.5in. / 67.3cm.	22in. / 55.9cm.	19in. / 48.3cm.	13in. / 33cm.
Large	Arm	D-302L	32in. / 81.3cm.	23in. / 58.4cm.	20in. / 50.8cm.	14in. / 35.6cm.

Lower Extremity		Model	Length	Circumference		
Size	Type	8 Chamber	Inseam	Ankle	Knee	Thigh/Groin**
Small	Full Leg	D-380S	23in. / 58.4cm.	21in. / 53.3cm.	25in. / 63.5cm.	30in. / 76.2cm.
Small Wide	Full Leg	D-380SW	26in. / 66cm.	24in. / 61cm.	29in. / 73.7cm.	35in. / 88.9cm.
Small Extra Wide	Full Leg	D-380SXW	24in. / 61cm.	28in. / 71.1cm.	35in. / 88.9cm.	40in. / 101.6cm.
Medium	Full Leg	D-380M	29in. / 73.7cm.	22in. / 55.9cm.	27in. / 68.6cm.	31in. / 78.8cm.
Medium Wide	Full Leg	D-380MW	29in. / 73.7cm.	25in. / 63.5cm.	31in. / 78.7cm.	35in. / 88.9cm.
Medium Extra Wide	Full Leg	D-380MXW	28in. / 71.1cm.	29in. / 73.7cm.	36in. / 91.4cm.	40in. / 101.6cm.
Large	Full Leg	D-380L	34in. / 86.4cm.	21in. / 53.3cm.	27in. / 68.6cm.	31in. / 78.7cm.
Large Wide	Full Leg	D-380LW	36in. / 91.4cm.	22in. / 55.9cm.	31in. / 78.7cm.	40in. / 101.6cm.

Upper Extremity		Model	Length*	Circumference		
Size	Type	8 Chamber	Axilla-Fingertips	Axilla	Elbow	Wrist
Medium	Arm	D-382M	26.5in. / 67.3cm.	22in. / 55.9cm.	19in. / 48.3cm.	13in. / 33cm.
Large	Arm	D-382L	32in. / 81.3cm.	23in. / 58.4cm.	20in. / 50.8cm.	14in. / 35.6cm.
Medium	Arm/Shoulder	D-383M	26.5in. / 67.3cm.	22in. / 55.9cm.	19in. / 48.3cm.	13in. / 33cm.
Large	Arm/Shoulder	D-383L	32in. / 81.3cm.	23in. / 58.4cm.	20in. / 50.8cm.	14in. / 35.6cm.

* Length corresponds to the entire sleeve length from axilla to fingertips ** The circumference of the upper rim of the sleeve